

**P315/PMC SURFACE DISINFECTANT 400 ml AMBRO-SOL**



**Features:**

Alcohol-based product designed for a quick and effective disinfection action of surfaces against bacteria, encapsulated viruses and yeasts. The product is colorless and versatile to use. Deeply cleans and disinfects any type of surface, evaporating quickly and leaving no residue.



**PRODUCT DATA**

VOLUME	400 ml
PACKAGING BOX	12 Pieces
BARCODE	8050038873208

**CHEMICAL-PHYSICAL PROPERTIES**


ASPECT	CAN CONTAINING UNDER PRESSURE LIQUID (AEROSOL)
COLOR	COLORLESS
ODOUR	CHARACTERISTIC
RELATIVE DENSITY	0,68 g/ml a 20°C
FLASH POINT	LOWER THAN 0 ° C - FLAMMABLE
FLASH POINT (LIQUID BASED)	< 23°C
PRESSURE A 20°C	5,5 bar ± 0,5 bar
INITIAL BOILING POINT (LIQUID BASED)	78°C

**LIQUID BASED ACTIVE PRINCIPLES**

**100 g of product contain:**

Isopropyl Alcohol	g	50,6
Ethyl Alcohol	g	5,6
Coformulants as required at	g	100,0

**HOW TO USE**



Shake the can well



direct the spray towards the surface to be disinfected, dispensing the product evenly at a distance of 15 cm. Wait for the complete evaporation of the product. It does not require rinsing except for surfaces that may come into contact with food.



Minimum contact time: 1 min.

**SURGICAL MEDICAL DEVICE (PMC)**

REG. n. 20981 of Italian Ministry of Health  
Authorization Owner: Ambro-Sol S.r.l.  
Via per Pavone del Mella, 21  
25020 Cigole (BS) - ITALY  
Tel: +39 030/9959674

**WARNINGS**



STORE IN A COOL AND DRY PLACE AT TEMPERATURES BETWEEN 5 °C AND 40 °C.

KEEP AWAY FROM HEAT SOURCES  
KEEP OUT OF THE REACH OF CHILDREN

**NOTES AND RECOMMENDATIONS**

IT'S RECOMMENDED TO VENTILATE THE ENVIROMENT BEFORE ACCESSING.

VALIDITY: 24 MONTHS FROM THE PRODUCTION DATE (CHECK AEROSOL BOTTOM)

AFTER USE, CLOSE THE CYLINDERS WITH THE CAP.

**BEFORE USE CONSULT THE SAFETY DATA SHEET OF THE PRODUCT**

This information, although believed to be reliable, must be considered indicative. Their use does not imply any liability on our part including the violation of any patents or licenses. It's the responsibility of the users to check in advance the suitability of the product for the specific use.